

Technical Data Sheet

EWB-ADR200, 200 ml FD

Order no.: 177 312

Product designation: EKASTU-Eye Wash Bottle ADR200, 200 ml FD

Applied standard: DIN EN 15154-4 (DIN = German Institute for Standardization)

Utilization: The first seconds after an accident, where the eyes are exposed to dangerous

substances, are decisive in keeping the injury to a minimum. Non plumbed-in eye wash units (eye wash bottles) are designed and intended to be kept in direct proximity to those people working in potentially dangerous areas. The principal duty

of these items of equipment is providing first aid.

Description: EKASTU Eye Wash Bottles are manufactured in accordance with DIN EN 15154-4

and are a Class I Medical Device, in accordance with the Directive 93/42 EEC.

Materials: The materials used are particularly kind of the skin and therefore have no negative

effect on human health.

Bottle: Low-density polyethylene (LDPE)
Cap: Polyethylene terephthalate (PET)
Connecting piece: Polyethylene-foam (PE-foam)
5 water channels strainer: Styrene-butadiene rubber (SBR)

Weight: approx. 245 g

Marking: Trade mark of manufacturer (LOGO)

Manufacturer/- address Product designation Applied standard

Consider the information brochure

LOT-No.

Expiry of shelf-life (MM/YYYY)

CE-Identification Storage conditions

Handling: Do not open the Eye Wash Bottle until immediately before use.

Storage conditions: Protect from sunlight, chemicals, heat and cold.

Observe storage conditions on the packaging. The eyewash bottle must be visually

inspected at regular intervals.

Warning reference: Please consider the information brochure!

User references: EKASTU Safety GmbH guarantees the indicated achievement according to class

Tel.:

Fax:

and type.

It is to be noted that the laboratory test values can considerably deviate from those, which are reached in practice. The user must read and understand all functional information. Use the Eye Wash Bottle with this information brochure, the relevant valid statutory regulations and the safety requirements of the profession associations

+49/7151/975099-0

+49/7151/975099-30

for example DGUV Information 204-022 'First Aid in Operation'.

info@ekastu.de

www.ekastu.de